

## Heritage MedCall Technical Specifications

### MODEL HM-GER-8 8' GERIATIRC CALL CORD

- Includes security clip
- Large pad for easy activation
- Safe for use in oxygenated environments
- 2-year warranty
- Latex-Free



#### **DESCRIPTION:**

The Heritage Model HM-GER-8 Geriatric Call Cord is a round, flattened momentary call cord pad that works great for residents who have limited hand dexterity. This call cord utilizes a mechanical switch in the pad that can be activated by slight pressure from the hand, arm or body, or can even be positioned under the user's chin.

#### **RESET:**

Emergency call station reset is accomplished by pressing the CHECK-IN or CANCEL button on the HM-S(D)PS panel. All emergency calls require this reset action to assure all calls are verified and all initiating locations are visited.

#### **INSTALLATION:**

Model HM-GER-8 is installed by fully inserting the 1/4" phone plug into the bed-cord wall plate. Removing the bed-cord from the wall plate will initiate an emergency call